

BIOBASE[®]

**UV/VIS SPECTROPHOTOMETER
BK-S360
User Manual**

BIOBASE GROUP

Version 2020.11

Content

Preface.....	2
1. Safety.....	2
2. About the Instrument.....	2
3. Symbols.....	2
I. Summary.....	3
1. Working Principle and Structure.....	3
2. Main Specifications.....	3
3. Main Functions.....	4
II. Install Instrument.....	5
1. Environment Required.....	5
2. Voltage Confirmation.....	5
3. Installation.....	5
III. Introduction of the Instrument.....	7
1. Introduction.....	7
2. Operation Panel.....	8
3. Keypad Description.....	8
IV. Instrument Operation.....	10
1. Software System.....	10
2. Basic Operation.....	11
3. Before Measurement.....	12
4. Measurement.....	12
5. Tools.....	24
V. Instrument Maintenance.....	26
1. Daily Maintain.....	26
2. Troubleshooting.....	26
3. Spare parts replacement.....	27

Preface

1. Safety

This instruments are designed according to the following safety documents: 《IEC 61010 — 1: 2001》 (Which is: 《Measurement, Control and Laboratory Electric Equipments Safety Requirement》 Part 1) and 《IEC 60601 — 1: 1988》 (Which is: 《Electric Apparatus of Medical Treatment》 Part 1).

2. About the Instrument

Scanning Series UV/Vis Spectrophotometer is an electrical test instruments which is widely used in the laboratories.

- Use Frequency: Intermittence
- Excessive Voltage (Current): No
- Pollution Class: Class 1

3. Symbols



Caution, Danger!



Caution, High Voltage!



Caution, Hot!



Ground



Fuse



Recall, this instrument will be called back by the appointed Electrical Treatment Department or by the original Manufacturer when wasted.

I. Summary

Scanning Series UV/Vis Spectrophotometer have the characters of wide range of wavelength, high sensitivity, powerful function, easy to use, simple structure and pretty figure. Besides these, the large LCD, High Precise A/D and easy to store RAM makes the instrument much more superior than other originals. It is widely used in Chemistry, Pharmaceuticals, Biochemical, metallurgy, Light Industry, Textile, Material, Environments, Medical, Education and some other fields. It is one of the most important instruments in Quality Control and an essential in normal laboratories.

1. Working Principle and Structure

Different matter has different but special absorbance wavelength point. Also, when at the fixed wavelength point, the absorbance has some relation to the substance's (Always transparent Solution) concentration and its thickness. The relation can be concluded as the following Formula which is called Lambert-Beer Law.

$$T = I/I_0$$

$$A = KCL = -\log I/I_0$$

A	Absorbance
C	Concentration of the Solution
K	Absorbance Coefficient of the Solution
L	The length of the Solution in the light path
I	The intensity of the light focused on the A/D after it permeate the solution to be measured.
I ₀	The intensity of the light focused on the A/D after it permeate the Solution.

Note: When test, the solvent is usually taken as the Reference Solution and its Transmittance is considered as 100%T. While the Transmittance of the sample to be tested is a relative value which is got comparing to that of the Reference.

2. Main Specifications

- Wavelength Range: 190-1100nm
- Stray Light: ≤0.05%T @ 220nm& 360nm
- Band Width: 1nm
- Photometric Range: 0-200%T,-0.3-3.0A
- Wavelength Accuracy: ±0.5nm
- WL. Repeatability: ≤0.2nm
- Photometric Accuracy: ±0.3%T
- Photometric Repeatability: ±0.2%T
- Stability: 0.001A/h @ 500nm
- Display: 480×272 Dots Matrix LCD
- Data Output: USB (2.0),USB (1.0)
- Dimension: 470×370×180
- Weight: 20kg

3. Main Functions

Scanning Series UV/Vis Spectrophotometer has functions below as main three parts.

- **Photometry**
The photometric mode can be switched between Absorbance, Transmittance and Energy. You can use single-point method to measure the concentration.
- **Quantitative**
Coefficient Method, Standard Curve Method and Input Method. The Regression Equation and the test result can be stored in the RAM and printed out.
- **Multiple Wavelength**
Test the Abs. and Trans of the same sample under different wavelengths, can test the data from 8 wavelengths.
- **Kinetics**
Test the sample changes in a period, get $\Delta A/t$, display the spectrum curve.
- **Spectrum**
Wavelength range 190-1100nm, set up the range, interval to test the max peak value curve of sample.
- **BIO**
Built-in 7 methods. Customers can also set the parameters. The data can be printed out.
- **Tools**
Users can set the favorable items to keep the instrument work in best conditions.

II. Install Instrument

1. Environment Required

To ensure the best performance, the following conditions are required:

- The best working temperature range is 16-35°C and the humidity is 45-80%.
- Keep it as far as possible away from the strong magnetic or electrical fields or any electrical device that may generate high-frequency fields.
- Set the unit up in an area that is free of dust, corrosive gases and strong vibrations.
- Remove any obstructions or materials that could hinder the flow of air under and around the instrument.
- The power requirement is 220±22V@50±1Hz or 110±11V@60±1 Hz.
- Use the appropriate power cord and plug into a grounded outlet.
- If the local voltage is not stable enough, a manostat is required.
- Be away from direct sunlight.

2. Voltage Confirmation



Be sure to set the instrument's voltage switch at your local power supply, or severe damage may occur!

3. Installation

Step 1: Check the packing list

Unpack the contents, check the materials with the Packing List. Any damage or Lost found, please contact us or the local dealer.

Step 2: Position

Place the instrument on the stable table carefully.

Step 3: Install printer (Optional)

Make sure the printer is power off. link the printer's data cable to the Instrument's parallel port.

Step 4: Link the power cord

Make sure the instrument's power switch is in the Off condition, link the power cord to the instrument and insert another end in the socket provided with a protective earth contact.

Step 5: Switch on the power

Check again. Make sure that all the links are right. Switch on the power. Then the instrument begins to self-test. After self-test and 15 minutes' pre-warm, it can work. The self-test includes the following steps:

Filter→Light Source Positioning→DarkCurrent→WLSelf-check→WEnergy→D Energy→ Room→ Battery Voltage(Fig.2-1)

After pre-warm, the instrument will ask the user to re-calibrate the system. (Fig. 2-2) Users can decide if they need to re-calibrate the system or not.




Fig. 2-1



Fig. 2-2

III. Introduction of the Instrument

1. Introduction

Appearance of the instrument (Fig. 3-1):



Fig. 3-1

- 1 — Lid of the room
- 2 — Operation Panel
- 3—Cell Holder
- 4— Rod
- 5— Fan
- 6— Fan Cover
- 7— Power Socket
- 8— Con-Tem Port
- 9— Power Switch
- 10—USB
- 11—Print Port
- 12 —USB Print Port

2. Operation Panel

Operation panel of Scanning UV/Vis Spectrophotometer(Fig.3-2):




Fig. 3-2

- 1 —LCD Displayer
- 2 —Keypad

3. Keypad Description

0 OPEN	9 WXYZ	Input Parameters, Wavelength and set 8 cell auto cell holder.
CLEAR		Delete the input value or stored data.
RETURN		Cancel operation or return to former menu
ZERO		Calibrate 100%T/0Abs
0 OPEN		Load the files saved in the RAM
SAVE		Save the test results.
START STOP		Start/Cancel Test

GOTO λ

Set Wavelength

PRINT

Print Test Result

ENTER

Confirm the input value or setting



Scroll the menu or data



Seek wavelength peak or set coordinate.

IV. Instrument Operation

1. Software System

Software System of Instrument (Fig. 4-1)




Fig. 4-1

2. Basic Operation

(1) Select Test Method

In the main menu, press numeric key to enter this mode directly.

(2) Set Wavelength

Press to set wavelength, use numeric key to input the value and press to confirm. When the wavelength changed to the point you set, the system calibrates 100%T/0Abs automatically.

(3) Set Parameters

In different menus, press to go into different parameter setting interface, input the values by numeric keys, press to confirm, press to return. Press to go into test interface and save settings.

(4) Set auto-cell holder(Optional Accessory)

Go into the auto-cell holder (In the setting menu) to make corresponding cell position at the light path. Then read the test data.

(5) Delete the Input Value

When input value, press to delete all the characters.

(6) Delete the test result and stored data

Press 、 to select in the test interface, press to delete the test result or stored data.

Press to delete all data if nothing selected.

(7) Calibrate 100%T/0Abs

Put the Reference in the light path, press to calibrate 100%T/0Abs.

(8) Measure Samples

In the test interface, put the samples to be tested in the light path, press to measure.

(9) Print the Test Result

In the test interface, press to print the test result.

(10) Store the Test Result

In the test interface, press to remind the store position, press again to input the file name and press to confirm.

(11) Load the Stored File

In the test interface, press  to go into file selecting interface, press 、 to choose the file you want, press  to open.



In different menu, users can only load the files with the corresponding postfix. List is the corresponding postfixes:

- Quantitative: *.qua
- Standard Curve: *.fit
- Wavelength Scan: *.wav
- Kinetics: *.kin
- DNA/Protein test: *.dna
- Multi-wavelength: *.mul

3. Before Measurement

(1) Self-check

Remove all the blocks in the light path and close the lid of the compartment. Switch on the power supply to begin the self-test.

(2) Pre-warming

After self-test, the instrument goes into pre-warm state. For accurate test, at least 30 minutes of warm up is required.

(3) Check the Cuvettes

The cuvettes must be clear and there’s no remains of the samples on the surface of it. Only Silicon (Quartz) cuvettes are permitted to be used in the range of less 300nm.

4. Measurement

(1) Basic Mode

Step 1, Go into test interface of Basic mode

In the main menu,press the key of  to choose “BASIC”followed with pressed to set up wavelength.




Fig. 4-2

Step 2, Set Test Mode

Press to set up (Fig. 4-2). Press to go into selecting mode, Press , to choose “TRANS”, “ABS” or “ENERGY” mode, and press to confirm. Press to set coefficient K (-9999.9-9999.9) to confirm. Press to confirm all settings.

Step 3, Set Wavelength

Press to set wavelength, input the value by the numeric keypad followed with pressed to confirm. Press to go into photometry data sheet.

Step 4, Calibrate 100%T/0Abs

Put the Reference in the light path and press to calibrate 100%T/0Abs.

Step 5, Sample Measurement

Put the sample to be measured in the light path, the result will be displayed on the screen automatically. Press to go into photometry data sheet. Repeat Step 4, 5, press to calculate the data into the data sheet, test all samples.

Step 6, Save/Open/Print Data

Press to save the result. Press , to choose store position, press to input file name, to confirm and save. Save photometry parameter and tested data. Fig. 4-3

Press , , to choose stored data, press to open and test directly.

Press to print data in data sheet.

BASIC DATA				500.0 nm -0.000 A			BASIC DATASAVE		
No.	DATA	K*DATA	PARAM	NAME	TIME	NOTICE			
1	0.000	0.0000	ABS WL:500.0 nm K :1.0000	SAMPLE	13-2-12 14:00	↑ ↓ ← →			
							SEL POS		
							SAVE		
							NAME		
							CLEAR		
							DEL CUR		
							1/1		

Fig. 4-3

(2) Quantitative

Step 1, Go into Quantitative Test Interface

In the main menu, press to choose “QTY”. The instrument display the last tested parameter.

QTY SAMPLE				500.0 nm			QTY SETUP		
No.	A	Conc.	PARAM	1. MODE	RECORD	ZERO			
1	0.411	0.3724	RECORD	2. FITTING	LINEAR	LINEAR			
2	0.605	0.5024	LINEAR	3. WL MODE	DOUBLE				
			UNIT: none	4. UNIT	none				
			WL1: 220.0	5. ROOM	MANUAL				
			SAMPLE: 2	6. M	1				
				7. N	2				

QTY SETUP		
1. WL1	220.0	RANGE
2. WL2	275.0	2
		1
		8
4. SAMPLE	2	2

Fig. 4-4

Step 2, Set Parameters

Press to enter the quantitative setting interface (Fig. 4-4). Press numeric key to select. Press to confirm. Press to enter standard sample test, input the standard sample concentration, press to confirm. After finished setting, press to test successively or press to input standard sample Abs. value.

P.S.: Coefficient method without displaying curve.

Step 3, Set up Standard Curve

Press to go into QTY Curveinterfaceto establish or select standard curve (Fig. 4-5), 3 methods are under your choice.




Fig. 4-5

Set up Standard Curve (In the setup menu):

Method 1: Standard Sample read (RECORD)

- 1) **Quantitative Method.** Press ,  ,  to choose standard sample read (RECORD).
 to confirm.
- 2) **Set Parameter.** According to the menu options, press numeric key to set up fitting method, wavelength method, units and room, press to confirm, then press to enter wavelength, sample no., to confirm.
- 3) **Sample Test.** Press to go into quantitative standard sample sheet.
- 4) **Set Standard Sample.** According to the notice, input the standard sample concentration (Fig. 4-7), to confirm, press  ,  to select the value. After finishing inputting concentration, press to enter Abs. test.
- 5) **Calibrate 100%T/0Abs.** Put the Reference Sample in the light path and press to calibrate 100%T/0Abs.
- 6) **Mark Standard Sample.** Put the sample corresponding with the concentration in the light path, press to read. Abs. The value will be displayed on the data sheet, according to notice to test successively.
- 7) **Sample Test.** Press to enter curve, press again to enter quantitative data sheet. Put the reference in the light path, then press . Put the sample in the light path, then press to read value.

Method 2: Use Standard Samples (INPUT)

- 1) **Quantitative Method.** Press  ,  to choose Standard Sample Input (INPUT), then press to confirm.
- 2) **Set Parameter.** According to the menu options, press numeric key to set up fitting method, wavelength method, units and room, press to confirm, then press to enter wavelength, sample no., to confirm.
- 3) **Sample Test.** Press to enter quantitative standard sample sheet.

4) **Set Standard Sample.** According to the notice, input the standard sample concentration (Fig. 4-7), to confirm, press ,  to select the value. After finishing inputting concentration, press to enter Abs. INPUT, the method is same as inputting concentration. Then press to enter curve picture, press again to enter quantitative data sheet to start test.

5) **Calibrate 100%T/0Abs.** Put the Reference Sample in the light path and press to calibrate 100%T/0Abs.

6) **Test Sample.** Put the test sample in the light path, press to read value.

Method 3: Coefficient method (Input Regression Equation)

1) **Quantitative Method.** Press , ,  to select coefficient method (PARAM). to confirm.

2) **Set Parameter.** According to the menu options, press numeric key to set up fitting method, wavelength method, units and room, press to confirm, then press to enter, input wavelength, photometry coefficient K1 and intercept K0 (Linear Fit) to confirm.

3) **Sample Test.** Press to enter quantitative data sheet. Fig. 4-6

QTY SETUP		QTY SETUP	
1. MODE	PARAM	1. WL1	500.0
2. FITTING	LINEAR		RANGE
3. WL MODE	DOUBLE		-9999.9
4. UNIT	NONE	4. K0	0.097
5. ROOM	MANUAL	5. K1	0.68
			+9999.9
			0.68

QTY SAMPLE		500.0 nm	
No.	A	Conc.	PARAM
1	0.411	0.3724	RECORD
2	0.605	0.5024	LINEAR
			UNIT:none
			WL1:220.0
			SAMPLE:2

Fig. 4-6

QTY SAMPLE				QTY SAMPLE			
No.	Conc.	A	NOTICE	NO.	C	A	PARAM
1	10		ZERO	1	10.000	0.104	RECORD
2	20		BLANK	2	20.000	0.210	LINEAR
3	30			3	30.000	0.305	UNIT:none
4	40		ENTER RECORD	4	40.000	0.408	WL1:220.0 SAMPLE

Fig. 4-7

Load the Stored Curves:

In the “Quantitative” interface, press to go into files selection interface. Use ,  to select the curve you need and press to load.

Users can press to check the curve.(Fig.4-5),press to return.

Step 4,Go into the sample measurement interface

In the “Quantitative Curve” interface, Press to go into the Quantitative Test interface(Fig. 4-8).

QTY SAMPLE			500.0 nm -0.001 A	
No.	A	Conc.	PARAM	
1	0.411	0.3724	WL: 500.0nm	
2	0.605	0.5024	K0: 0.097	
			K1: 0.67	
			A=K1*C+K0	

Fig. 4-8

Step 5, Calibrate 100%T/0Abs.

Put the Reference Sample in the light path, press to calibrate 100%T/0Abs.

Step 6, Measure Samples

Place the sample to be tested in the light path, press to measure. Then the test result will display in the data sheet. Repeat this step to finish measuring all the samples.(Fig. 4-9)

QTY SAMPLE			500.0 nm -0.001 A		QTY SAMPLE			500.0 nm	
No.	A	Conc.	PARAM		No.	A	Conc.	PARAM	
1			WL: 500.0nm		1	0.411	0.3724	RECORD	
2			K0: 0.097		2	0.605	0.5024	LINEAR	
			K1: 0.67					UNIT:none	
			A=K1*C+K0					WL1:500.0	
								SAMPLE:2	

Fig. 4-9

Step 7, Print Data

Press  to print the result.

Step 8, Delete Data

Press ,  to move the cursor to select the value you don't need, then press  to delete the value.

Step 9, Store the data

After measurement, press  to store the data. Input the file name by the numeric keypad and press  to store.

Double Wavelength: $\Delta A = M \cdot A_1 - N \cdot A_2$

Triple Wavelength: $\Delta A = A_2 - A_4$

$$A_4 = \frac{(\lambda_1 - \lambda_2)A_3 + (\lambda_2 - \lambda_3)A_1}{\lambda_1 - \lambda_3}$$

(3) Kinetics

Step 1, Go into Kinetics Test Interface

In the main menu, press  to select "KINS" (Fig. 4-12).




Fig. 4-12

Step 2, Setup Parameters

Press function key  to set up parameters, input wavelength, test time, delay, interval, record range and room.  to confirm.

Press numeric key  to choose test mode Abs. or Trans.

Press numeric key  to input test wavelength, range 190nm-1100nm.

Press numeric key  to input test time, the unit is second.

Press numeric key  to input delay time, the delay time cannot be calculated to data sheet.

Press numeric key

5
JKL

 to input interval, the data in the data sheet read according to the interval.

Press numeric key

6
MNO

 to input the sample Abs. range, which is -4-4.

Press numeric key

7
PQRS

 to choose room, standard is manual, default is auto.

After finished, press

START
STOP

 to enter spectrum interface.

Step 3, Calibrate 100%T/0Abs

Put the Reference Sample in the light path, press

ZERO

 to calibrate 100%T/0Abs.

Step 4, Measure the sample

Put the sample to be tested in the light path and press

START
STOP

 to begin the test.

Step5, Data Process

After test, press

GOTO λ

 to enter data process interface. (Fig.4-13)




Fig. 4-13

Press

1

 to recover original data.

Press

2
ABC

 to change coordinate, test time and range.

Press

3
DEF

 to enter data sheet, display test data A and dA/dt. Fig. 4-14

Press

◀

,

▶

 to display the point value successively.

KINS DATA			
TIME	A	dA/dt	NOTICE
0	0.000		
10	0.000	0.0000	
20	0.000	0.0000	
30	-0.000	-0.0000	
40	0.000	0.0000	
			1/1

Fig. 4-14

Step 6, Print curve

Press to print curve and data sheet.

Step 7, Save and Open

After scan, press to save curve, parameters and data sheet, press numeric key to input file name, to save. to open stored file. , to choose file, to open.

(4) Multi-Wavelength Test

Step 1, Go into multi-wavelength test interface

In the main menu, press numeric key to choose “MULTWL”. (Fig. 4-10)

MULTWL DATA				MULTWL SETUP		
No.	WL	A	NOTICE	1. MODE	TRANS	TRANS
1	275.0	0.002	ZERO	2. WL SUM	2	ABS
2	220.0	0.502	BLANK	3. INPUT WL	275.0 220.0	
			START RECORD			
			1/2			

Fig. 4-10

Step 2, Setup

Press , to set up test mode, press , to choose Trans. and Abs., to confirm.

Press numeric key to set up no. of wavelength, can set up 8 wavelengths at most, choose the wavelength no. under your choice.

Press numeric key to set up wavelength, then key to confirm.

After finished, press to start testing.

Step 3, Calibrate 100%T/0Abs

Press the sample in the light path, according to notice, press to calibrate 100%T/0Abs.

Step 4, Data Test

Place the Reference Sample in the light path, press to test sample, results will be displayed in the data sheet. (Fig. 4-11)

Step 5, Save Data

Press to save data, press , to choose store data, press again to input file name, to confirm. Press , , to choose stored file, to open and test directly. to print the tested data.

MULTWL DATA				BASIC DATASAVE			
No.	WL	A	NOTICE	NAME	TIME	NOTICE	
1	275.0	0.002	ZERO	SAMPLE	13-2-12 14:00	↑ ↓ ← →	
2	220.0	0.502	BLANK			SEL POS	
			START RECORD			SAVE NAME	
			1/2			CLEAR	
						DEL CUR	

Fig. 4-11

Step 6, Delete

Press to delete the current data.

(5) Spectrum

Step 1, Go into spectrum scan

In the main menu, press numeric key to choose "SPEC". (Fig. 4-15)

SPEC		650.0nm		SPEC SETUP		
T	100.0	PARAM		1. MODE	TRANS	TRANS
		MODE: TRANS		2. LIGHT	AUTO	ABS
0	230.0	LIGHT: AUTO		3. INTERVAL	0.1 nm	ENERGE
		INTERVAL: 0.1		4. XRANGE	650.0- 230.0	
	WL	SPEED: MID		5. YRANGE	0.0/100.0	
	650.0	GOTO		6. SPEED	LOW	
		PROCESS				

Fig. 4-15

Step 2, Set up scan parameters

Function key to set up parameters, according to notice, set up Mode, Source Light, Interval, Scan Range, Record Range and Scan Speed, to confirm, after set up, press to enter spectrum test interface.

Press numeric key , mode: Trans., Abs., Energy, press , key to choose test mode.

Press to confirm.

Press numeric key

2
ABC

, source light: Auto, W lamp, D lamp, press , to choose test source light. Press

ENTER

 to confirm.

Press numeric key

3
DEF

, set interval:0.1, 0.2, 0.5, 1, 2, 5nm, press , to choose scan interval. Press

ENTER

 to confirm.

Press

4
GHI

, input scan range: 190-1100nm, input the range under the practical meet. Press

ENTER

 to confirm.

Press

5
JKL

, input record range: Abs., A value is -4-4; Trans., T value is 0-200%; Energy, E value is 0-100%. Press

ENTER

 to confirm.

Press

6
MNO

, set scan speed: Fast, Middle, Slow, press , to choose scan speed, press

ENTER

 to confirm.

After finishing settings, press

START
STOP

 to enter spectrum test interface.

Step 3, Scan baseline

Put the reference in the light path,

ZERO

 to scan baseline,

START
STOP

 to cancel scan.

Step 4, Scan Sample

Put the sample in the light path,

START
STOP

 to scan sample (Fig. 4-16a),

START
STOP

 to cancel scan.




Fig. 4-16a




Fig. 4-16b

Step 5, Data Process

After scan, press GOTO λ to enter data process (Fig. 4-16b). Press numeric key 1, recover original spectrum; Press numeric key 2 ABC to change X,Y coordinate. Press numeric key 3 DEF to enter data sheet.

Step 6, Print Curve and Peak Value

Press PRINT to print curve and data sheet.

Step 7, Save Data

Press SAVE to choose store position, press SAVE again to save, input curve name and press ENTER to confirm.

Step 8, Open and Delete Curve

Press 0 OPEN, ▲, ▼ to choose the curve you want and ENTER to open, press CLEAR to delete stored curve.

(6) BIO

Step 1, Go into DNA/Protein Test

In the Main menu, press 6 MNO to choose “BIO”.(Fig. 4-17)




Fig. 4-17

DNA-1 DATA			DNA-1 SETUP		
NO.	WL	A	NOTICE		RANGE
	320.0		ZERO	1. BACKGROUND ON (320nm)	-9999.9
	260.0		BLANK	2. K1 49.100	
	230.0			3. K2 3.4800	+9999.9
				4. K3 183.00	
				5. K4 75.800	
Abs	RATIO		START		
DNA	CONC		RECORD		
PROTEIN	CONC				
			1/1		

Fig. 4-18

Step 2, Choose Method

Press numeric key 1 to choose DAN-1, Fig. 4-18

Step 3, Setup

Press to set coefficients, input all the values of 1-5 by numeric keypad according the indication and press to confirm. Press to enter test interface.

Step 4, Calibrate 100%T/0Abs

Place the Reference Sample in the light path, press to calibrate 100%T/0Abs.

Step 5, Sample Measure

Place the sample to be test in the light path, press to measure. The result will display in the data sheet.(Fig. 4-18).

Step 6, Delete Data

Press to delete the data.

Step 7, Store the Data

After measurement, press to store the result. Input the file name by the numeric keypad, press to confirm.

5. Tools

In the main menu,press to select “Tools”.(Fig. 4-19)

TOOLS			
1. W LAMP	ON	44H	W LAMP
2. D LAMP	ON	36H	<input type="button" value="ON"/>
3. DARKCURRENT			<input type="button" value="OFF"/>
4. WL CAL			
5. LANGUAGESENGLISH			
6. CRT		13-03-01 11:00	
7. SWPOINT		339.0nm	
8. REST			

Fig. 4-19

1. W Lamp

Press and use , to choose “W Lampon or off” and press to confirm.



Users can switch Off the W Lamp when test in the wavelength range of 190-339nm, so as to prolong the life time of W lamp.

2. D2 Lamp

Press  and use ,  to choose “D Lamp on or off” then press  to confirm.



Users can switch Off the D2 Lamp when test in the wavelength range of 340-1100nm, so as to prolong the life time of D2 lamp.

After some period of usage(About 1000h), the energy would drop off to a low level,and the test data in the UV range would get instable and deviate a lot from the experiential value. Then the D2 lamp should be replaced to a new one.After replacement, users should calibrate the wavelength again.

3. Dark Current

Press  to choose “DARKCURRENT” and press  to confirm. Then the system begins to refresh Dark Current.

Note: During the course, opening the lid of the compartment is prohibited.

4. Wavelength Calibration

Press  to re-search D lamp energy curve 656.1nm. The instrument will self-calibrate when switch on.

Note: During the course, opening the lid of the compartment is prohibited.

5. Languages Selection

Press  and use ,  to choose “Languages” and press  to confirm.

6. Clock Setup

Press  to choose “Clock Setup”, then press function key 1,2,3,4,5 to set up data (Fig. 4-20).

TOOLS				
1. W LAMP	ON	44H	1. Year:	13
2. D LAMP	ON	36H	2. Month:	03
3. DARKCURRENT			3. Day:	01
4. WL CAL			4. Hour:	11
5. LANGUAGES	ENGLISH		5. Minute:	00
6. RTC	13-03-01	11:00		
7. SWPOINT	339.0nm			
8. REST				

Fig. 4-20

7. Change Switch Point

Press numeric key  to choose “SWPOINT”. The default value is 339nm.

8. RESET

Press numeric key  to initial files, empty all ROM files. **NOTICE: Be Cautious to use this function.**

9. BASELINE

Press numeric key  to choose “BASELINE”. This function is for re-establish Baseline.

V. Instrument Maintenance

To keep the instrument work in good condition, constant maintain is needed.

1. Daily Maintain

(1) Check the Room

After measurement, the cuvettes with sample solutions should be taken out of the compartment in time. Or the volatilization of the solution would make the mirror go moldy. Users must pay more attention to the corrosive sample and liquid easy to volatilize. Any solution remains in the compartment should be wipe off immediately.

(2) Surface Clean

The cover of the instrument is with paint. Please use wet towel to wipe off the drips on the surface immediately. Organic solution is forbidden to be used to clean the cover. Please wipe off the dirt on the cover timely.

(3) Clean the Cuvettes

After every test or after a solution change, the cuvettes should be cleaned carefully, or the remains on the surface would cause measuring error.

2. Troubleshooting

(1) Dark Current Error when Self-check

Possible Cause

- Open the lid of the compartment during the course of self-test.

Solution

Close the lid of the compartment and switch on the power again.

(2) No Response After Power On

Possible Cause

- Bad contact in power supply
- Fuse melt

Solution

Improve the contact
Replace a new fuse

(3) Printer Not Work, Printer Error

Possible Cause

- No power supply
- Bad contact in power supply
- Bad contact of the data cable

Solution

Switch on the power supply
Improve the contact
Improve the contact

(4) No Stable Reading

Possible Cause

- No enough pre-warm
- Glass cuvettes used in UV Range
- No stable Sample
- Much higher sample concentration
- Low voltage or unstable power supply
- Lights defect
- Light used up

Solution

Increase the pre-warm time
Use Silicon Cuvettes.
Improve the sample
Dilute the sample
Improve the power condition
Replace a new lamp
Replace a new lamp

(5) Worse Repeatability

Possible Reason

- Unstable sample

Solution

Improve the sample

- Cuvettes polluted Clean the cuvettes

(6) Incorrect Reading

Possible Reason	Solution
· Dark Current Error	Re-get the dark current
· Worse matching of the cuvettes	Improve the matching of the cuvettes

3. Spare parts replacement

(1) Replace the Fuse



Danger! Be sure to switch off the power and unplug the socket before replacement!

Step 1, Tools Preparation

Prepare a 3×75 Flat Blade screwdriver.

Step, Switch Off the Power Supply

Switch off the power supply, and unplug the socket.

Step 3, Take Out the Fuse Seat

Take out the fuse seat by the screwdriver(Fig.5-1).




Fig. 5-1

Step 4, Replace a New Fuse

Pick out the spare fuse and replace it to the working position(Fig.5-2).




Fig. 5-2

Step 5, Reset the Fuse Seat

Replace the fuse seat in the power socket.

Step 6, Switch On the Power

Plug the socket and switch on the power.

(2) Replace Lamps



High temperature ! Wait 20 minutes before open the lamp chamber after power off to avoid scald!

Step 1, Tools Preparation

Prepare a 6×150mm Cross Blade screwdriver and a pair of glove.

Step 2, Power Off

Switch off the power supply and unplug the socket.

Step 3, Open the Cover

Unscrew the 4 screws indicated in Fig.5-3(Each side with 2 screws)and remove the cover.




Fig. 5-3

Step 4, Open the Cover of the Light Chamber

Unscrew the 2 screws on the light chamber cover and remove it(Fig.5-4).




Fig. 5-4

Step 5, Replace the D2 Lamp

Unscrew the 2 screws on the D2 Flange (No.1 in Fig. 5-5), unplug the connector in the power board (No. 2) and remove the D2 lamp. Draw on the cotton glove and replace a new lamp. Fix the 2 screws and plug the connector again.




Fig. 5-5

Step 6, Replace W Lamp



Remember the direction of the filament before pull out the W lamp. Be sure that the new lamp's filament is in the same direction as before.

Pull out the defected W lamp and draw on the cotton glove. Insert the new W lamp as deep as possible on the lamp seat. Be sure to keep the filament in the same direction as the old one face.




Fig. 5-6

Step 7, Adjust the Position of the W Lamp

Switch on the power,(the Switch Mirror should be placed to the position as Fig. 5-7 indicates). Observe the entrance facula,and it should in the center of the entrance hole(Fig.5-7). If the facula deviate to Left or Right, then loosen the No.1 screws in Fig. 5-8 and move the lamp seat to Left or Right until it focus on the center of the slot. Then fix the screws. If the facula deviate to Up and Down, then loosen the No.2 screws in Fig. 5-8 and move the lamp seat Up and Down until the facula focus on the center of the slot. Then fix the No. 2 screws again.




Fig. 5-7




Fig. 5-8

Step 8, Finish

Reset the cover of the light chamber and fix the screws. Reset the cover of the instrument and fix the screws. Recover the Pole in the compartment, then the course finished.

(3) Replace the Battery



Be sure to switch off the power supply and unplug the socket before open the Bottom Cover!

Step 1, Prepare the Tools

Prepare a 6×150mm Cross Blade Screwdriver.

Step 2, Switch Off the Power Supply

Switch off the power supply and unplug the socket.

Step 3, Open the Bottom Cover Plate

Unscrew the 13 screws indicated in Fig.5-9 then remove the bottom plate.



Fig. 5-9

Step 4, Replace the Battery

Pick out the old battery and replace a new one(Fig. 5-10).




Fig. 5-10

Step 5, Finish

Recover the bottom plate and fix the 13 screws, then the course finishes.

Appendix 1: Consumable Sheet

Deuterium Lamp

Tungsten Lamp

10mm Cuvette Glass (4PCS/BOX)

10mm Cuvette Quartz (4PCS/BOX)

BIOBASE GROUP

2# building, No.9 Gangxing Road, High-tech Zone, Jinan City, Shandong Province,
China

Tel: +86-531-81219803/01

Fax: +86-531-81219804

Inquiry: export@biobase.com

Complaints: customer_support@biobase.cc

After-sales service: service_sd@biobase.cc; service_ivd@biobase.cc

Web: www.biobase.cc/www.meihuatrade.com / www.biobase.com