

BIOLGLX®

Manufacturer and Supplier of Laboratory Products

www.BiologixGroup.com

BIOLGLX® V. 2021




USER MANUAL

Centrifuge Tube Racks

Series 90 - XXX

- *Place the rack in a safe and stable place, preferably with enough space to avoid knocks with nearby elements during use.*
- *Insert the different sample tubes or test tubes into the available holes in a vertical position, always making sure it remains firm and completely vertical.*



Maintenance and Cleaning

Microcentrifuges, PP centrifuges, and Cardboard Centrifuge Made of high strength plastic, chemically resistant to alcohols and mild organic solvents **DO NOT REQUIRE MAINTENANCE.**

For cleaning

- Use a cloth moistened with 70% ethanol and go over the Surface.

Applicable to:

| Cat # | Well | Color | Sterile | Unit | Total |
|---------|---------|-----------------|---------|----------------------|-----------|
| 90-8009 | 80-well | Assorted Colors | No | 5 /Bag, 4 Bags/Case | 20 /Case |
| 90-1536 | 36-well | White | No | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-5016 | 16-well | White | No | 5 /Bag, 20 Bags/Case | 100 /Case |
| 90-1550 | 50-well | Blue | No | 1 /Bag, 20 Bags/Case | 20 /Case |
| 90-5025 | 25-well | Blue | No | 1 /Bag, 20 Bags/Case | 20 /Case |