

BIOLGLX®

Manufacturer and Supplier of Laboratory Products

www.BiologixGroup.com

BIOLGLX® V. 2021




USER MANUAL

Centrifuge Tube Racks

Series 90 - XXX

- *Place the rack in a safe and stable place, preferably with enough space to avoid knocks with nearby elements during use.*
- *Insert the different sample tubes or test tubes into the available holes in a vertical position, always making sure it remains firm and completely vertical.*



Maintenance and Cleaning

Microcentrifuges, PP centrifuges, and Cardboard Centrifuge Made of high strength plastic, chemically resistant to alcohols and mild organic solvents **DO NOT REQUIRE MAINTENANCE.**

For cleaning

- Use a cloth moistened with 70% ethanol and go over the Surface.

Applicable to:

Cat #	Well	Color	Sterile	Unit	Total
90-8009	80-well	Assorted Colors	No	5 /Bag, 4 Bags/Case	20 /Case
90-1536	36-well	White	No	5 /Bag, 20 Bags/Case	100 /Case
90-5016	16-well	White	No	5 /Bag, 20 Bags/Case	100 /Case
90-1550	50-well	Blue	No	1 /Bag, 20 Bags/Case	20 /Case
90-5025	25-well	Blue	No	1 /Bag, 20 Bags/Case	20 /Case